

Peel-Away Introducer
14 cm Sheath
Introducer Kit
7 F – 16 F


Product Highlights

- Proprietary materials improve insertion characteristics and reduce vessel trauma
- Close tolerance extrusion and proprietary tipping process improves tracking on a guidewire
- Reliable peeling characteristics
- Sheath can be totally occluded without kinking to prevent air inspiration
- Di-Lock™ feature secures dilator in sheath during insertion

Ordering Information

Contents: Peel-Away Sheath, Di-Lock™ Dilator, 12 cc syringe, 18 ga. XTW Needle, and 50 cm Guidewire with 3 mm “J”

Model Number	Size (F)	Maximum Guidewire Diameter (in)	Sheath Usable Length (cm)	Units per Box
405154	7 F-TW	,038	14	1
405145	8 F-TW	,038	14	1
405146	8 F	,038	14	1
405147	9 F	,038	14	1
405149	10 F	,038	14	1
405104	6 F	,038	14	5
405108	7 F	,038	14	5
405112	8 F	,038	14	5
405129	8 F-TW	,038	14	5
405116	9 F	,038	14	5
405118	9,5 F	,038	14	5
405120	10 F	,038	14	5
405122	10,5 F	,038	14	5
405124	11 F	,038	14	5
405128	12 F	,038	14	5
405132	13 F	,038	14	5
405144	16 F	,038	14	5



Customer Support: 46-8-474-4756

Brief Summary: Prior to using these devices, please review the Instructions for Use for a complete listing of indications, contraindications, warnings, precautions, potential adverse events and directions for use. Devices depicted may not be available in all countries. Check with your St. Jude Medical representative for product availability in your country.

Unless otherwise noted, ™ indicates that the name is a trademark of, or licensed to, St. Jude Medical or one of its subsidiaries. ST. JUDE MEDICAL and the nine-squares symbol are registered and unregistered trademarks and service marks of St. Jude Medical, Inc. and its related companies. ©2015 St. Jude Medical, Inc. All rights reserved.

SJM-CAG-1115-0055 | Item approved for international use only.